

# PHANTOM BACK



The Phantom Back™ is designed to target active wheelchair users, who require back support and need maximum freedom for upper body movement.


The Phantom Back™ fits rigid wheelchair frames without push handles and replaces standard back support upholstery. It offers maximum seat-to-back infinite angle adjustment within the range, providing real individual fit. The well designed back support panel and cushion promote lateral and posterior postural support in pelvic region, needed for fatigue reduction, thus enhancing upper body function.

The accessibility of hardware interface, the flexibility to adjust and it's lightweight make the Phantom Back™ exceptional in meeting the needs of an active user.

Phantom Series		STYLE: Phantom Back
Backs		
Stock Sizes – Same Day Shipment		
ADP CODE: SEBFM0070		
Part Number	Size W x H	
PB127	12" x 7"	
PB137	13" x 7"	
PB147	14" x 7"	
PB157	15" x 7"	
PB167	16" x 7"	



Cover, Foam Layers and Aluminum Hardware



## Product Highlights:

### ✓ Height Adjustments

Infinite adjustment along wheelchair posts. Hardware and slotted holes in the back panel provide 3" adjustment range increments of 1", to overcome cross brace obstacle on the wheelchair posts.

### ✓ Angle Adjustments

In relation to the wheelchair back support posts, the hardware accommodates up to 15 degrees anterior and up to 15 degrees of posterior infinite angle adjustment within this range.

### ✓ Wheelchair Posts Interface

Standard clamp hardware will properly fit 1" tubing diameter. For wheelchairs with 7/8" or 3/4" tubing diameter, insert spacers are included in the hardware.

### ✓ Space for Storage

Pocket provided on back for storage of personal items.

### ✓ Width Adjustments

Hardware allows width adjustment range of 3/4". For example 16" wide back would fit a wheelchair frame from 15-3/4" to 16-1/2".

### ✓ Seat Depth Adjustments

Optimal seat depth of 1-1/4" behind wheelchair posts available with standard posts. Maximum seat depth is 2-1/2" behind the wheelchair posts, when anterior back angle is 15 degrees; minimum seat depth is 0" behind the wheelchair posts, when posterior angle is set to 15 degrees. For shorter seat depth than 0", a custom spacer is available to set the seat depth up to 2" in front of wheelchair posts.

### ✓ Product Weight: 2.5 lbs (16"W x 7"H)

**User Weight:** 250 lbs

**Warranty:** 3 years limited on back, 90 days on cover, and 1 year on hardware

Lateral Dimensions - Phantom		
PHANTOM, Heights 14, 16, 18, 20		
Chair Width	DIM B (")	DIM A (")
12"	9.5	1.4
13"	10.5	1.4
14"	11.5	1.4
15"	12.5	1.4
16"	13.5	1.4
17"	14.5	1.4
18"	15.5	1.4

