

LUNAR BACK

LUNAR



The Lunar Back™ was designed for users who require aggressive lateral support and more complex positioning. Lateral support options such as fixed or adjustable swing aways make this back the most versatile and truly modular on the market. Lunar's sleek and attractive tailored look will appeal to a wide variety of users!


Product Highlights:

- ✓ Lightweight, sturdy aluminum shell
- ✓ Attractive cover design
- ✓ Replaceable foam
- ✓ Available in 3 back style heights (low back, mid-back, and full back)
- ✓ Available in 3 lateral styles (no lateral, basic lateral, and ultra lateral)
- ✓ Optional lateral supports, available in 4 styles: Fixed, Height Adjustable, Fixed Offset and Adjustable Swing Away
- ✓ Pins adjustable to accommodate up to 2" difference in seat width
- ✓ **Product Weight:** 5.9 lbs (18"W x 18"H Basic)
Capacity: 250 lbs
Warranty: 3 years limited on back and 90 days on cover, 1 year on hardware



Lateral Styles (Non, Basic, Ultra)





Lateral Dimensions - LUNAR					
NO LATERAL, Heights 11-26"			BASIC, Height 16"		
Chair Width	DIM B (")	DIM A (")	Chair Width	DIM B (")	
14"	11.75	2.64	14"	12.58	4.32
16"	13.75	2.64	16"	14.58	4.32
18"	15.75	2.64	18"	16.58	4.32
20"	17.75	2.64	20"	18.58	4.32
22"	19.75	2.64	22"	20.58	4.32
ULTRA, Heights 10-26"					
Chair Width	DIM B (")	DIM A (")			
14"	13.75	6.34			
16"	15.75	6.34			
18"	17.75	6.34			
20"	19.75	6.34			
22"	21.75	6.34			



- ✓ Fully adjustable for height, depth and angle of recline
- ✓ Hardware fits 3/4", 7/8" and 1" tubing and can be installed within a few minutes with the right tools
- ✓ Quick release lower latch levers allows the backrest to be easily removed.