

# LUNAR HD BACK




The Lunar Back Heavy Duty™ was designed for users weighing up to 500lbs who require aggressive lateral support and more complex positioning. Reinforcement ribs and extra-large shell allow this back to maintain its function in any high load condition. Lunar's sleek and attractive tailored look will appeal to a wide variety of users!

**Lunar Heavy Duty** **STYLE: Mid Support**

Stock Sizes – Same Day Shipment

ADP CODE: SEBFM3055

Part Number	Size W x H	Description	Part Number	Size W x H	Description	Part Number	Size W x H	Description
MBSHDBL2216	22" x 16"	Basic	MBSHDBL2616	26" x 16"	Basic	MBSHDBL3016	30" x 16"	Basic
MBSHDUL2216	22" x 16"	Ultra	MBSHDUL2616	26" x 16"	Ultra	MBSHDUL3016	30" x 16"	Ultra
MBSHDBL2218	22" x 18"	Basic	MBSHDBL2618	26" x 18"	Basic	MBSHDBL3018	30" x 18"	Basic
MBSHDUL2218	22" x 18"	Ultra	MBSHDUL2618	26" x 18"	Ultra	MBSHDUL3018	30" x 18"	Ultra
MBSHDBL2416	24" x 16"	Basic	MBSHDBL2816	28" x 16"	Basic	MBSHDBL3216	32" x 16"	Basic
MBSHDUL2416	24" x 16"	Ultra	MBSHDUL2816	28" x 16"	Ultra	MBSHDUL3216	32" x 16"	Ultra
MBSHDBL2418	24" x 18"	Basic	MBSHDBL2818	28" x 18"	Basic	MBSHDBL3218	32" x 18"	Basic
MBSHDUL2418	24" x 18"	Ultra	MBSHDUL2818	28" x 18"	Ultra	MBSHDUL 3218	32" x 18"	Ultra

**Lunar Heavy Duty** **STYLE: Full Support**

Stock Sizes – Same Day Shipment

ADP CODE: SEBFM3060

Part Number	Size W x H	Description	Part Number	Size W x H	Description	Part Number	Size W x H	Description
FBSHDBL2220	22" x 20"	Basic	FBSHDBL2620	26" x 20"	Basic	FBSHDBL3020	30" x 20"	Basic
FBSHDUL2220	22" x 20"	Ultra	FBSHDUL2620	26" x 20"	Ultra	FBSHDUL3020	30" x 20"	Ultra
FBSHDBL2222	22" x 22"	Basic	FBSHDBL2622	26" x 22"	Basic	FBSHDBL3022	30" x 22"	Basic
FBSHDUL2222	22" x 22"	Ultra	FBSHDUL2622	26" x 22"	Ultra	FBSHDUL3022	30" x 22"	Ultra
FBSHDBL2224	22" x 24"	Basic	FBSHDBL2624	26" x 24"	Basic	FBSHDBL3024	30" x 24"	Basic
FBSHDUL2224	22" x 24"	Ultra	FBSHDUL2624	26" x 24"	Ultra	FBSHDUL3024	30" x 24"	Ultra
FBSHDBL2420	24" x 20"	Basic	FBSHDBL2820	28" x 20"	Basic	FBSHDBL3220	32" x 20"	Basic
FBSHDUL2420	24" x 20"	Ultra	FBSHDUL2820	28" x 20"	Ultra	FBSHDUL3220	32" x 20"	Ultra
FBSHDBL2422	24" x 22"	Basic	FBSHDBL2822	28" x 22"	Basic	FBSHDBL3222	32" x 22"	Basic
FBSHDUL 2422	24" x 22"	Ultra	FBSHDUL2822	28" x 22"	Ultra	FBSHDUL3222	32" x 22"	Ultra
FBSHDBL2424	24" x 24"	Basic	FBSHDBL2824	28" x 24"	Basic	FBSHDBL3224	32" x 24"	Basic
FBSHDUL2424	24" x 24"	Ultra	FBSHDUL2824	28" x 24"	Ultra	FBSHDUL3224	32" x 24"	Ultra

### Product Highlights:

- ✓ Lightweight, sturdy aluminum shell
- ✓ Reinforcement ribs
- ✓ Attractive cover design
- ✓ Replaceable foam
- ✓ Available in 2 lateral styles (basic lateral and ultra lateral)
- ✓ Available in 2 back style heights (mid-back and full back)
- ✓ Optional lateral supports, available in 4 styles: Fixed, Height
- ✓ Adjustable, Fixed Offset and Adjustable Swing Away
- ✓ Pins adjustable to accommodate up to 2" difference in seat width
- ✓ **Product Weight:** - 8.55 lbs (22"W x 18"H Basic)
- Capacity:** 500 lbs
- Warranty:** 3 years limited on back and 90 days on cover, 1 year on hardware



Lunar HD showing Inner Cover



- ✓ Fully adjustable for height, depth and angle of recline
- ✓ Hardware fits 3/4", 7/8" and 1" tubing and can be installed within a few minutes with the right tools
- ✓ Quick release lower latch levers allows the backrest to be easily removed.

