

Orion III



TILTING WHEELCHAIR
ORION III

 **FUTURE**
MOBILITY HEALTHCARE INC.



A. DUAL CROSS GAS CYLINDERS

DUAL CROSS GAS CYLINDERS LOCATED UNDERNEATH THE CHAIR, ALLOWS FOR A GENTLE, SLOW TILT FOR THE CLIENT, LOCKABLE AT INFINITE ANGLES

B. TILT IN SPACE

ORION III OFFERS A TRUE 45 DEGREE OF TILT FOR OPTIMAL POSITIONING FOR THE INDIVIDUAL WITH AN OPTIONAL 30 DEGREE OF RECLINE

C. ATTENDANT FOOT TILT

ATTENDANT FOOT TILT ALLOWS THE WHEELCHAIR ATTENDANT TO EASILY TILT WHEELCHAIR BY PRESSING ON THE PEDAL TO TILT UP OR DOWNWARDS.

D. POWER TILT

POWER TILT IS AVAILABLE FOR EITHER THE ATTENDANT OR USER TO TILT THE WHEELCHAIR WITHOUT ANY ADDED MANUAL FORCE.

E. ATTENDANT TRIGGER TILT

ALL BASE MODELS OF ORION III ARE EQUIPPED WITH THE TRIGGER TILT SYSTEM. THIS ALLOWS ATTENDANTS OF THE USER TO TILT THE WHEELCHAIR EASILY.

F. ADJUSTABLE AXLE AND SEAT HEIGHT

REAR AXLE MOUNT CAN BE ADJUSTED FOR GREATER STABILITY AND MOVEMENT



WEIGHT SHIFT CENTER OF MASS
TILT TRACE PATH GIVING A MORE
STABLE FINAL TILTED POSITION



FEATURES ON THE ORION III WHEELCHAIR

- TILT OPTION
- RECLINE OPTION
- LOW SEAT-TO-FLOOR HEIGHT
- STANDARD SEATING SYSTEM CONSISTING OF: TWO WAY CONTOURED PADDED FIXED HEADREST, ORION BACK WITH 1" DEPTH ADJUSTABILITY AND PRISM IDEAL CUSHION
- ALUMINUM SEAT PAN
- MULTI ADJUSTABLE AXLE PLATE
- REMOVABLE REAR WHEELS WITH STAINLESS STEEL OR PLASTIC COATED HANDRIMS
- ADJUSTABLE ARM HEIGHTS, GEL PADDING
- 70 DEGREE PIN STYLE FRONT RIGGING
- ALUMINUM OR COMPOSITE FOOTPLATES
- QUICK FLIP-UP ANTI TIPPERS
- LOCKING COLLAR
- HEEL LOOPS

CUSTOM ANYTHING!

CUSTOMER COMES FIRST. FUTURE MOBILITY HEALTHCARE IS A CANADIAN MANUFACTURER OF WHEELCHAIRS, BACKS, CUSHIONS AND ACCESSORIES WITH OUR PARTS MANUFACTURED IN HOUSE. THIS GIVES US THE ABILITY TO CUSTOMIZE ANYTHING NEEDED FOR SPECIFIC CUSTOMER REQUESTS.














THE ORION III IS THE MOST INNOVATIVE AND AFFORDABLE TILT-IN-SPACE WHEELCHAIR ON THE MARKET TODAY. DEVELOPED FOR THOSE INDIVIDUALS REQUIRING ADJUSTABILITY TO REDISTRIBUTE THEIR BODY WEIGHT FOR OPTIMAL POSITIONING PURPOSES, THE ORION III OFFERS 3-45 DEGREE OF TILT, OPTIONAL -9, +33 DEGREE OF INFINITE RECLINE, STANDARD ON-CHAIR ADJUSTABLE SEATING SYSTEM, INCREDIBLE RANGE OF SEAT-TO-FLOOR HEIGHTS AND CREATIVE STREAMLINED DESIGN.

Orion III

WWW.FUTUREMOBILITY.CA

ORION III

SPECIFICATIONS:

 SEAT WIDTH 16"- 24"	 SEAT TO FLOOR 13"- 20"	 REAR WHEEL 12", 16", 20", 22" OR 24" COMPOSITE URETHANE
 SEAT DEPTH 15"- 24"	 OVERALL WIDTH NARROW SEAT WIDTH + 9.5" REGULAR SEAT WIDTH + 10.0	 CASTER SIZE 4", 5", 6", 8" URETHANE
 FRONT RIGGING SWINGAWAY DUAL SWINGAWAY OR DUAL ELEVATING LEGRESTS	 ADJUSTABLE BACK ANGLE -7°, 0°, 7°, 14°, 21°	 BACK STYLE LOW PUSH CANES W/STROLLER BAR HIGH PUSH CANES W/STROLLER BAR 16", 18", 20" ORION BACK
 ARM STYLE ADJ. HEIGHT T-STYLE ADJ. HEIGHT PIN RELEASE T-STYLE ADJ. HEIGHT FLIP BACK	 PRODUCT WEIGHT 65 LBS (29 KG)	 WEIGHT CAPACITY 300 LBS (159 KG)
		 ADP Ontario CODE WM0001195

WARRANTY:

THE ORION III IS WARRANTED TO BE FREE FROM DEFECTS IN MATERIALS AND WORKMANSHIP FOR A PERIOD OF TWO (2) YEARS FROM THE DATE OF PURCHASE AND LIFETIME FOR THE FRAME. UPHOLSTERY FOR 90 DAYS, UNDER NORMAL USE BY ORIGINAL PURCHASER.

POWER TILT BATTERIES ARE WARRANTED FOR SIX (6) MONTHS WHILE THE CONTROL BOX, BATTERY CHARGER, ACTUATOR AND HAND CONTROL ARE WARRANTED FOR ONE (1) YEAR.



FUTURE MOBILITY HEALTHCARE BELIEVES IN CONTINUED ENHANCEMENT, RESEARCH, INVENTION AND IMPROVEMENT OF MOBILITY PRODUCTS. OUR PRODUCTS CONTINUE TO IMPROVE TO MEET THE NEEDS OF OUR CUSTOMERS FOR A BETTER QUALITY OF LIFE.

FUTURE MOBILITY HEALTHCARE INC. 3223 ORLANDO DRIVE, MISSISSAUGA, ONTARIO, L4V 1C5

1-888-737-4011 FAX: 1-866-670-3377

WWW.FUTUREMOBILITY.CA

1 BUFFALO RIVER PLACE BUFFALO, NEW YORK, 14210, USA

1-855-9-FUTURE FAX: 716-783-9236,

WWW.FUTUREMOBILITY.COM

ALL RIGHTS RESERVED

PRINTED IN CANADA.

INNOVATIVE SEATING AND MOBILITY SOLUTIONS



EUROPEAN CONFORMANCE

