

PRISM SUPER BACK FOR KIDS



The Prism Super Back™ for Kids is a junior size seat back which provides children excellent pressure relief. Extras, such as the positioning option, the ease of attachment and the unique method of adjusting for a wide range of angles, make this back the most suitable for special pediatric needs.


Prism Series Backs		STYLE: Super Back for Kids	
Stock Sizes - Same Day Shipment			
ADP CODE: SEBFM0080			
Part Number	Size W x H	Part Number	Size W x H
KSB1012	10" x 12"	KSB1414	14" x 14"
KSB1212	12" x 12"	KSB1616	16" x 16"
Non-Stock Sizes – Shipment within 3 business days			
ADP CODE: SEBFM0085			
Part Number	Size W x H	Part Number	Size W x H
KSB1112	11" x 12"	KSB1412	14" x 12"
KSB1114	11" x 14"	KSB1416	14" x 16"
KSB1214	12" x 14"	KSB1512	15" x 12"
KSB1216	12" x 16"	KSB1514	15" x 14"
KSB1312	13" x 12"	KSB1612	16" x 12"
KSB1316	13" x 16"	KSB1614	16" x 14"
Custom Sizes - Shipment within 7 business days			
ADP CODE: SEBFM0090			
Part Number	Size W x H		
KSBCUS	Custom		



Back view of Foam Layer and Aluminum Shell

Product Highlights:

- ✓ Lightweight, sturdy aluminum 5052 grade shell
- ✓ Full 5" range of forward and rear adjustment (angle adjustment forward 15" and back 20")
- ✓ Pins adjustable to accommodate up to 2" difference in seat width (standard, 0.5" and 1" longer)
- ✓ Attractive cover design
- ✓ Cane clamps hardware included
- ✓ Replaceable foam
- ✓ Lumbar pad and contour pad are standard
- ✓ Quick push button release for easy installation
- ✓ **Product Weight:** 3.2 lbs (12"W x 14"H)
- ✓ **User Weight:** 150 lbs
- ✓ **Warranty:** 3 years limited on back, 90 days on cover, and 1 year on hardware



Lateral Dimensions - Kids Back		
KIDS BACK, Heights 10, 12, 13, 14, 16		
Chair Width	DIM B (")	DIM A (")
12"	9.6	4.5
13"	10.6	4.5
14"	11.6	4.5
15"	12.6	4.5
16"	13.6	4.5



- ✓ Fully adjustable for height, depth and angle of recline
- ✓ Hardware fits 3/4", 7/8" and 1" tubing and can be installed within a few minutes with the right tools
- ✓ Quick release lower latch levers allows the backrest to be easily removed.