


# PRISM SUPER CUSHION FOR KIDS



The Super Cushion for Kids™ is a Junior size cushion created with three layers of polyurethane, high density and memory foam, which provides children with excellent pressure relief. The anatomically designed firm foam base with its aggressive leg trough and abductor contour, provides superior pelvic stability and leg positioning. A central core of pressure relieving visco-elastic insert under the IT's (Ischial Tuberosity), allows for excellent immersion for bony prominences.



#### CUSHION THICKNESS

- A- Height at adductor = 3 3/4"
- B- Height at abductor = 4"
- C- Height at rear IT = 2 1/2"
- D- Leg Trough = 3 1/2"
- E- Height at rear trochanter = 3 1/2"

#### Prism Series Cushions STYLE: Super Cushion for Kids Stock Sizes – Same Day Shipment ADP CODE: SESFM0090

Part Number	Size W x L	Part Number	Size W x L
KSC1010	10" x 10"	KSC1416	14" x 16"
KSC1212	12" x 12"	KSC1614	16" x 14"
KSC1414	14" x 14"	KSC1616	16" x 16"

#### Non-Stock Sizes – Shipment within 3 business days ADP CODE: SESFM0095

Part Number	Size W x L	Part Number	Size W x L
KSC1012	10" x 12"	KSC1216	12" x 16"
KSC1014	10" x 14"	KSC1410	14" x 10"
KSC1016	10" x 16"	KSC1412	14" x 12"
KSC1210	12" x 10"	KSC1610	16" x 10"
KSC1214	12" x 14"	KSC1612	16" x 12"

#### Custom Sizes - Shipment within 7 business days ADP CODE: SESFM0100


Part Number	Size W x L
KSCCUS	Custom




Cushion Layers and Inner Cover

### Product Highlights:

- ✓ High density layered foam for a consistent, long lasting product
- ✓ Standard water resistant outer cover, made of black Meshtex urethane
- ✓ Waterfall front edge for foot propelling comfort
- ✓ Radius cut to match shape of back upholstery
- ✓ Pressure relief
- ✓ Optional 4-way stretch, water resistant inner liner
- ✓ 4-way stretch, breathable outer cover
- ✓ **Product Weight:** 2.5 lbs (14"W x 14"L)
- ✓ **User Weight:** 150 lbs
- ✓ **Warranty:** Two (2) years limited on cushion and 90 days on the cover.



Pressure Reading Scale in mmHG



Pressure Map	Supreme Gel
Weight	150 Lbs/68 Kg
Duration	1 hr
Average Pressure	24 mmHg

