

Adjustable Drop Base

Future Mobility Healthcare manufactures and markets products that increase the freedom and mobility to allow for a better quality of life. Renown for quality and service we bring innovation to all of our products for the home healthcare market.

The Prism Adjustable Drop Base was designed for clients to provide a stable base of support on the seat cushion and adjust in width in 1/2" increments, up to a 2" maximum adjustment.

The Prism Adjustable Drop Base lowers the seat cushion height for an improved fit under tables or desks and also allows for easier foot propulsion. The Prism Adjustable Drop Base is made of aluminum and is lightweight. Its unique design allows the cushion to breathe easier, yet still provide the stability to the client.



PROUDLY CANADIAN



Adjustable Drop Base



SPOTLIGHT ON

- Flexible sizes
- Attractive design

WARRANTY

One (1) years limited.

Stock sizes—same day shipments

16-18 X 15	ADB16-1815
16-18 X 16	ADB16-1816
16-18 X 17	ADB16-1817
16-18 X 18	ADB16-1818
16-18 X 19	ADB16-1819
16-18 X 20	ADB16-1820

Non-stock sizes—shipments within 3 business days

14-16 X 16	ADB14-1616
14-16 X 18	ADB14-1618
19-21 x 16	ADB19-2116
19-21 X 17	ADB19-2117
19-21 X 18	ADB19-2118
19-21 X 19	ADB19-2119
19-21 X 20	ADB19-2120

Custom sizes—shipments within 7 business days

ADBCUS

PRODUCT WEIGHT

6 lbs without cushion

USER WEIGHT

250 lbs

Future M  bility Healthcare

3223 Orlando Drive, Mississauga, ON, L4V 1C5, 905-671-1661, 888-737-4011, fax 866-670-3377 www.future-mobility.com
Future Mobility Healthcare Inc. Orion and Prism are trademarks of Future Mobility Healthcare Inc.