

Serene Prime



HCPCS: E2601/D2 General Use Cushion

THE PRIME CUSHION was designed for users that require mild postural support and offers gentle contours, helping to promote correct seating posture. A seat cushion made of medium density, high resilient, polyurethane foam, that is lightweight, comfortable and with good pressure management.

PRODUCT HIGHLIGHTS:

- ✓ Standard water proof cover
- ✓ Radiused front edge for foot propelling
- ✓ Light weight cushion
- ✓ Non slip bottom on washable, incontinent cover
- ✓ Promotes mild postural support
- ✓ 4 way stretch, water resistant inner liner
- ✓ Medium density, high resilient, polyurethane foam




Non slip Bottom on washable, incontinent cover




Polyurethane foam, anatomically shaped

Product Weight: 1.5 Lbs / 0.68 Kg (16"W x20"L)



Pressure Reading Scale in mmHG



User Weight 250 Lbs / 113 Kg
 Duration 1 hr
 Average Pressure 36 mmHg

SERENE CUSHION		STYLE: PRIME	
Stock Sizes – Same Day Shipment			
Part #	Size W x L	Part #	Size W x L
PR1616	16" x 16"	PR1818	18" x 18"
PR1618	16" x 18"	PR2016	20" x 16"
PR1816	18" x 16"	PR2018	20" x 18"

SERENE CUSHION		STYLE: PRIME	
Non-Stock Sizes – Shipment within 3 business days			
Part #	Size W x L	Part #	Size W x L
PR1620	16" x 20"	PR2216	22" x 16"
PR1820	18" x 20"	PR2218	22" x 18"
PR2020	20" x 20"	PR2220	22" x 20"
		PR2416	24" x 16"
		PR2418	24" x 18"
		PR2420	24" x 20"