



**FUTURE**  
 MOBILITY PRODUCTS INC.  
 WWW.FUTUREMOBILITY.COM  
 TOLL FREE 1-855-938-8873



1 Buffalo Rive Place Buffalo, NY, 14210

<b>Account Information</b>		Order Date: _____
Account #: _____	Business Name: _____	P.O.#: _____
Ship to Address: _____	Fax: _____	Phone: _____
Contact: _____		City: _____
		State: _____
		Zip Code: _____

**The HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions. Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy. Prices are subject to change without notice**

Item	Parts #	E-Code	Retail Price
<b>1. FRAME TYPE</b>			
<input type="checkbox"/> Galaxy Lite Folding Chair	103	K0005	\$2,195.00
<b>2. FRAME FINISH</b>			
<input type="checkbox"/> Wet Black (Standard)	103C01	-	No Charge
<input type="checkbox"/> Silver	103C02	-	No Charge
<input type="checkbox"/> Blue	103C03	-	No Charge
<input type="checkbox"/> Green	103C04	-	No Charge
<input type="checkbox"/> Orange	103C05	-	No Charge
<input type="checkbox"/> Purple	103C06	-	No Charge
<input type="checkbox"/> Red	103C07	-	No Charge
<input type="checkbox"/> Yellow	103C08	-	No Charge
<b>3. SEAT WIDTH</b>			
<input type="checkbox"/> 12" Seat Width	103SW12	-	\$75.00
<input type="checkbox"/> 13" Seat Width	103SW13	-	No Charge
<input type="checkbox"/> 14" Seat Width	103SW14	-	No Charge
<input type="checkbox"/> 15" Seat Width	103SW15	-	No Charge
<input type="checkbox"/> 16" Seat Width	103SW16	-	No Charge
<input type="checkbox"/> 17" Seat Width	103SW17	-	No Charge
<input type="checkbox"/> 18" Seat Width	103SW18	-	No Charge
<input type="checkbox"/> 19" Seat Width	103SW19	-	No Charge
<input type="checkbox"/> 20" Seat Width	103SW20	E2201	\$425.00
<b>4. SEAT DEPTH</b>			
<input type="checkbox"/> 12" Seat Depth	103SD12	-	\$75.00
<input type="checkbox"/> 13" Seat Depth	103SD13	-	No Charge
<input type="checkbox"/> 14" Seat Depth	103SD14	-	No Charge
<input type="checkbox"/> 15" Seat Depth	103SD15	-	No Charge
<input type="checkbox"/> 16" Seat Depth	103SD16	-	No Charge
<input type="checkbox"/> 17" Seat Depth	103SD17	-	No Charge
<input type="checkbox"/> 18" Seat Depth	103SD18	-	No Charge
<input type="checkbox"/> 19" Seat Depth	103SD19	-	No Charge
<input type="checkbox"/> 20" Seat Depth	103SD20	E2203	\$425.00

**Standard Chair Includes:**  
 Lightweight Folding Frame  
 Elliptical Heavyduty Crossbrace Tubing  
 Seat Width (12", 14"-20")  
 Seat Depth (12", 15"-20")  
 Seat-To-Floor Height Front (13" - 19-1/2")  
 Seat-To-Floor Height Rear (12-1/2" -18-1/4")  
 Rear Wheel Options (20", 22" or 24")  
 Handrim Options (Stainless Steel or Plastic Coated)  
 Caster Options (4",5",6" or 7")  
 Slip on Seat Upholstery  
 Slip on Back Upholstery  
 60 Degree Swingaway Front Riggings  
 Angle Adjustable Push Canes

Colors Option in Stellato Wheelchair



**FOOTNOTES:**

- A. Seat angle is set at 3 degrees. See service manual for details
- I. Only available as a 12"x12" chair with STF 13.5" to 16" on a hemi frame. Not available with 22" rear wheels. Must choose Superchair footplates.
- K. Allow a few extra days for powder coating time for custom colors

Item	Parts #	E-Code	Retail Price
<b>5. SEAT TO FLOOR HEIGHT</b>			
<input type="checkbox"/> 13.0" Seat to Floor (20" Wheels & 4" Casters) Low frame	103STF45	K0056	\$100.00
<input type="checkbox"/> 14.0" Seat to Floor (20" Wheels & 4" Casters) Low frame	103STF46	K0056	\$100.00
<input type="checkbox"/> 14.5" Seat to Floor (20" Wheels & 4" Casters) Low frame	103STF47	K0056	\$100.00
<input type="checkbox"/> 15.0" Seat to Floor (20" Wheels & 4" Casters) Low frame	103STF48	K0056	\$100.00
<input type="checkbox"/> 15.0" Seat to Floor (20" Wheels & 5" Casters) Low frame	103STF03	K0056	\$100.00
<input type="checkbox"/> 15.5" Seat to Floor (20" Wheels & 5" Casters) Low frame	103STF04	K0056	\$100.00
<input type="checkbox"/> 15.5" Seat to Floor (22" Wheels & 4" Casters) Low frame	103STF49	K0056	\$100.00
<input type="checkbox"/> 16.0" Seat to Floor (20" Wheels & 4" Casters) Low frame	103STF50	K0056	\$100.00
<input type="checkbox"/> 16.0" Seat to Floor (22" Wheels & 5" Casters) Low frame	103STF05	K0056	\$100.00
<input type="checkbox"/> 16.5" Seat to Floor (20" Wheels & 5" Casters) Low frame	103STF51	K0056	\$100.00
<input type="checkbox"/> 17.0" Seat to Floor (20" Wheels & 6" Casters) STD frame	103STF52	-	No Charge
<input type="checkbox"/> 17.0" Seat to Floor (22" Wheels & 5" Casters) STD frame	103STF53	-	No Charge
<input type="checkbox"/> 17.5" Seat to Floor (20" Wheels & 6" Casters) STD frame	103STF54	-	No Charge
<input type="checkbox"/> 17.5" Seat to Floor (22" Wheels & 5" Casters) STD frame	103STF55	-	No Charge
<input type="checkbox"/> 18.0" Seat to Floor (20" Wheels & 6" Casters) STD frame	103STF56	-	No Charge
<input type="checkbox"/> 18.0" Seat to Floor (22" Wheels & 6" Casters) STD frame	103STF57	-	No Charge
<input type="checkbox"/> 18.0" Seat to Floor (24" Wheels & 6" Casters) STD frame	103STF20	-	No Charge
<input type="checkbox"/> 18.5" Seat to Floor (22" Wheels & 6" Casters) STD frame	103STF58	-	No Charge
<input type="checkbox"/> 18.5" Seat to Floor (24" Wheels & 6" Casters) STD frame	103STF59	-	No Charge
<input type="checkbox"/> 19.0" Seat to Floor (22" Wheels & 7" Casters) STD frame	103STF60	-	No Charge
<input type="checkbox"/> 19.5" Seat to Floor (22" Wheels & 7" Casters) STD frame	103STF61	-	No Charge
<b>6. REAR WHEEL</b>			
<input type="checkbox"/> Mag Wheels	103RW09	-	No Charge
<input type="checkbox"/> Spoke Wheels	103RW10	-	\$255.00
<b>7. REAR WHEEL TIRE</b>			
<input type="checkbox"/> Airless (Solid Rubber)	103RT1	-	No Charge
<input type="checkbox"/> Pneumatic	103RT2	E2211	\$60.00
<b>8. CHAMBER OPTIONS</b>			
<input type="checkbox"/> 0° Camber	103CB1	-	No Charge
<input type="checkbox"/> 3° Camber	103CB2	-	No Charge
<b>9. HANDRIM (Available with 20", 22" or 24" Rear Wheels)</b>			
<input type="checkbox"/> Chrome	103HR1	-	No Charge
<input type="checkbox"/> Plastic Coated	103HR2	K0059	\$90.00
<input type="checkbox"/> No Handrims	103HR3	-	No Charge
<b>10. ARM TYPE</b>			
<input type="checkbox"/> T-Style Adjustable Height Pin Release - Youth	103AR05	E0973	\$210.00
<input type="checkbox"/> T-Style Adjustable Height Pin Release - Adult	103AR06	E0973	\$210.00
<input type="checkbox"/> T-Style Fixed Height	103AR02	E0973	\$210.00
<input type="checkbox"/> Full Length Adj. Height Flipback	103AR09	E0973	\$215.00
<input type="checkbox"/> Full Length Fixed Height Flipback	103AR10	E0973	\$195.00
<input type="checkbox"/> Swingaway Padded Tubular Arm	103AR01	E0973	\$175.00
<b>11. ARM PAD TYPE</b>			
<input type="checkbox"/> Desk length Vinyl 10" (pair)	103AP1	-	No Charge
<input type="checkbox"/> Full length Vinyl 14" (pair)	103AP2	-	No Charge
<input type="checkbox"/> Contour Gel 17" (pair)	103AP3	-	No Charge
<b>12. BACK POST</b>			
<input type="checkbox"/> Angle Adj 15 deg Bend w/Push Handle (low)	103PC3	-	No Charge
<input type="checkbox"/> Angle Adj 15 deg Bend w/Push Handle (high)	103PC4	-	\$90.00

SEAT TO FLOOR MATRIX										
GALAXY LITE FRONT STF HEIGHT AVAILABILITY										
REAR WHEEL	4" CASTER									
20"	13"	14"	14.5"	15"	15.5"	16"				
22"										
REAR WHEEL	5" CASTER									
20"				15"	15.5"	16.5"				
22"						15"	17"	17.5"		
REAR WHEEL	6" CASTER									
20"							17"	17.5"	18"	
22"									18"	18.5"
24"									18"	18.5"
REAR WHEEL	7" CASTER									
22"										19"

**FOOTNOTES:**

B. Can not use a 8" caster with this wheelchair

C. Short fork is used to achieve these front seat to floor heights

D. For front seat to floor of 14" and above the long fork is used



FULL LENGTH FIXED HEIGHT FLIP BACK



FULL LENGTH HEIGHT ADJ. FLIP BACK



GEL FULL LENGTH HEIGHT ADJ. FLIP BACK

Item	Parts #	E-Code	Retail Price
------	---------	--------	--------------

**13. BACK POST HEIGHT**

<input type="checkbox"/> 16" Post Height	NA	-	No Charge
<input type="checkbox"/> 17" Post Height	NA	-	No Charge
<input type="checkbox"/> 18" Post Height	NA	-	No Charge
<input type="checkbox"/> 19" Post Height	NA	-	No Charge
<input type="checkbox"/> 20" Post Height	NA	-	No Charge
<input type="checkbox"/> 21" Post Height	NA	-	No Charge
<input type="checkbox"/> 22" Post Height	NA	-	No Charge
<input type="checkbox"/> 23" Post Height	NA	-	No Charge
<input type="checkbox"/> 24" Post Height	NA	-	No Charge

**14. BACK OPTIONS**

<input type="checkbox"/> No Back	101BK49	-	Standard
<input type="checkbox"/> Slip-on Padded Back Upholstery	103BK48	-	No Charge
<input type="checkbox"/> Adjustable Tension Back (High 16-20")	TBH16-20	E2611	\$295.00
<input type="checkbox"/> Adjustable Tension Back (Low 13-17")	TBL16-20	E2611	\$295.00
<input type="checkbox"/> Prism Multiback High (16")	MBH14-20	E2611	\$349.00
<input type="checkbox"/> Prism Multiback Low (13")	MBL14-20	E2611	\$330.00
<input type="checkbox"/> Prism Basic Back 14-20" wide	BB14-20	E2620	\$615.00
<input type="checkbox"/> Prism Basic Back 22-26" wide	BB22-26	E2621	\$750.00
<input type="checkbox"/> Prism Basic Back 27-32" wide	BBCUS	E2621	\$775.00
<input type="checkbox"/> Prism Ultra Back 14-20" wide	UB14-20	E2620	\$615.00
<input type="checkbox"/> Prism Ultra Back 22-26" wide	UB22-26	E2621	\$750.00
<input type="checkbox"/> Prism Ultra Back 27-32" wide	UBBCUS	E2621	\$775.00
<input type="checkbox"/> Prism Air Back 14-20" wide	AB14-20	E2615	\$655.00
<input type="checkbox"/> Prism Air Back 22-26" wide	AB22-26	E2616	\$755.00
<input type="checkbox"/> Prism Air Back 27-32" wide	ABCUS	E2616	\$815.00
<input type="checkbox"/> Prism TrueFitt Back 14-20" wide	TFB14-20	E2615	\$625.00
<input type="checkbox"/> Prism TrueFitt Back 22-26" wide	TFB22-26	E2616	\$725.00
<input type="checkbox"/> Prism TrueFitt Back 27-32" wide	TFBCUS	E2616	\$775.00
<input type="checkbox"/> Prism Max Ultra Back 14-20" wide	MUB14-20	E2620	\$645.00
<input type="checkbox"/> Prism Max Ultra Back 22-26" wide	MUB22-26	E2620	\$695.00
<input type="checkbox"/> Prism Max Ultra Back 27-32" wide	MUBCUS	E2620/21	\$725.00
<input type="checkbox"/> Prism Max Air Back 14-20" wide	MAB14-20	E2620	\$725.00
<input type="checkbox"/> Prism Max Air Back 22-26" wide	MAB12-32	E2620	\$755.00
<input type="checkbox"/> Prism Max Air Back 27-32" wide	MABCUS	E2620/21	\$775.00
<input type="checkbox"/> Prism Max TrueFitt Back 14-20" wide	MTFB14-20	E2620	\$675.00
<input type="checkbox"/> Prism Max TrueFitt Back 22-26" wide	MTFB22-32	E2620	\$725.00
<input type="checkbox"/> Prism Max TrueFitt Back 27-32" wide	MTFBCUS	E2620	\$775.00
<input type="checkbox"/> Lunar Back Low No Lateral on chair 14-20"	LBSNL14-20	E2620	\$595.00
<input type="checkbox"/> Lunar Back Low Basic Lateral on chair 14-20"	LBSBL14-20	E2620	\$670.00
<input type="checkbox"/> Lunar Back Mid No Lateral on chair 14-20"	MBSNL14-20	E2620	\$645.00
<input type="checkbox"/> Lunar Back Mid Basic Lateral on chair 14-20"	MBSBL14-20	E2620	\$720.00
<input type="checkbox"/> Lunar Back Mid Ultra Lateral on chair 14-20"	MBSUL14-20	E2620	\$795.00
<input type="checkbox"/> Lunar Back Full No Lateral on chair 14-20"	FBSNL16-22	E2620/21	\$695.00
<input type="checkbox"/> Lunar Back Full Basic Lateral on chair 14-20"	FBSBL16-22	E2620/21	\$770.00
<input type="checkbox"/> Lunar Back Full Ultra Lateral on chair 14-20"	FBSUL16-22	E2620/21	\$845.00
<input type="checkbox"/> Lunar Back Custom on chair 14-20"	LUNARCUS	E2620/21	\$ Call



Prism Max Ultra Back with 6" laterals



Prism Ultra Back with 6" laterals



Prism Basic Back with 3" laterals

**FOOTNOTES:**

F. Sizes 14,16,18,20W x 10H or 12H

G. Sizes 14,16,18,20W x 14H,16H or 18H

H. Sizes 16,18,20,22W x 20H,22H or 24H

(refer to lunar backs order form for


specific sizes)

**15. BACK PINS (Prism Backs only)**

<input type="checkbox"/> Longer Pins for Back H/Ware 1" Growth	103BP1	-	\$50.00
<input type="checkbox"/> Longer Pins for Back Hardware 2" Growth	103BP2	-	\$50.00

Item	Parts #	E-Code	Retail Price
------	---------	--------	--------------

16. SEAT CUSHION			
<input type="checkbox"/> No Cushion	101SC20	-	Standard
<input type="checkbox"/> Black Nylon Seat Upholstery	103SC19	-	Standard
<input type="checkbox"/> Serene Prime Cushion 16-20" wide	PR16-24	E2601	\$100.00
<input type="checkbox"/> Serene Paramount Cushions 16-20" wide	PA16-20	E2605	\$300.00
<input type="checkbox"/> Prism Ideal Cushion 16-20" wide	IC16-20	E2601	\$100.00
<input type="checkbox"/> Prism Ideal Cushion 22-26" wide	IC22-26	E2601	\$175.00
<input type="checkbox"/> Prism Ideal Cushion 27-32" wide	ICCU	E2601/02	\$225.00
<input type="checkbox"/> Prism Elite Cushion 16-20" wide	EC16-20	E2603	\$225.00
<input type="checkbox"/> Prism Elite Cushion 22-26" wide	EC22-26	E2603	\$300.00
<input type="checkbox"/> Prism Elite Cushion 27-32" wide	ECCUS	E2603/04	\$375.00
<input type="checkbox"/> Prism Supreme Cushion 16-20" wide	SC16-20	E2607	\$365.00
<input type="checkbox"/> Prism Supreme Cushion 22-26" wide	SC22-26	E2607	\$495.00
<input type="checkbox"/> Prism Supreme Cushion 27-32" wide	SCCUS	E2607/08	\$595.00
<input type="checkbox"/> Prism Supreme II Plus 16-20" wide	SCP16-20	E2607	\$385.00
<input type="checkbox"/> Prism Supreme II Plus 22-26" wide	SCP22-26	E2607	\$425.00
<input type="checkbox"/> Prism Supreme II Plus 27-32" wide	SCPCUS	E2607/08	\$475.00
<input type="checkbox"/> Prism Supreme Air Cushion 16-20" wide	SAC16-20	E2624	\$475.00
<input type="checkbox"/> Prism Supreme Air Cushion 22-26" wide	SAC22-26	E2625	\$550.00
<input type="checkbox"/> Prism Supreme Air Cushion 27-32" wide	SACCUS	E2624/25	\$595.00
<input type="checkbox"/> Prism Supreme Gel Cushion 16-20" wide	SGC16-20	E2607	\$395.00
<input type="checkbox"/> Prism Supreme Gel Cushion 22-26" wide	SGC22-26	E2607	\$425.00
<input type="checkbox"/> Prism Supreme Gel Cushion 27-32" wide	SGCUS	E2607/08	\$475.00



Supreme II Plus



17. FRONT RIGGING			
-------------------	--	--	--

<input type="checkbox"/> 70° Pin Style	103FR01	-	No Charge
<input type="checkbox"/> 80° Pin Style	103FR02	-	No Charge
<input type="checkbox"/> 60° Swingaway	103FR03	-	No Charge
<input type="checkbox"/> 70° Swingaway	103FR04	-	No Charge
<input type="checkbox"/> 90° Swingaway (Not available w/ 8" casters)	103FR05	-	No Charge
<input type="checkbox"/> 60° Dual Rotation	103FR06	-	\$230.00
<input type="checkbox"/> 70° Dual Rotation	103FR07	-	\$230.00
<input type="checkbox"/> 80° Dual Rotation	103FR08	-	\$230.00
<input type="checkbox"/> 90° Dual Rotation	103FR09	-	\$230.00
<input type="checkbox"/> Swingaway Elevating Legrest (pair)	103FR10	E0990	\$275.00
<input type="checkbox"/> Dual Swingaway Elevating Legrest (pair)	103FR25	E0990	\$290.00
<input type="checkbox"/> No Front Rigging	103FR11	-	No Charge



18. FOOTPLATE			
---------------	--	--	--

<input type="checkbox"/> Composite Flip Up Footplates Small	103FP01	-	No Charge
<input type="checkbox"/> Composite Flip Up Footplates Medium	103FP02	-	No Charge
<input type="checkbox"/> Composite Flip Up Footplates Large	103FP03	-	No Charge
<input type="checkbox"/> Aluminum Flip Up Footplates Small	103FP05	-	No Charge
<input type="checkbox"/> Aluminum Flip Up Footplates Medium	103FP06	-	No Charge
<input type="checkbox"/> Aluminum Flip Up Footplates Large	103FP07	-	No Charge
<input type="checkbox"/> Aluminum Angle Adj Flip Up Footplates Small	103FP10	K0040	\$120.00
<input type="checkbox"/> Aluminum Angle Adj Flip Up Footplates Medium	103FP11	K0040	\$120.00
<input type="checkbox"/> Aluminum Angle Adj Flip Up Footplates Large	103FP12	K0040	\$120.00
<input type="checkbox"/> One Piece Footboard	103FP19	-	\$160.00



19. FOOTREST/LEGREST ACCESSORIES			
----------------------------------	--	--	--

<input type="checkbox"/> Padded Calf Strap	103FA03	K0038	\$45.00
<input type="checkbox"/> Foam Pad on Footplate (pair)	103FA07	-	\$45.00
<input type="checkbox"/> Foam Pad on One Piece Footboard	103FA06	-	\$50.00
<input type="checkbox"/> Heel Loops	103FA01	K0951	\$40.00
<input type="checkbox"/> No Heel Loops	103FA02	-	No Charge

Item	Parts #	E-Code	Retail Price
------	---------	--------	--------------

20. WHEEL LOCKS			
<input type="checkbox"/> Push to Lock	103WL1	-	No Charge
<input type="checkbox"/> Push to Lock + Wheel Lock Extension	103WL3	E0961	\$60.00
<input type="checkbox"/> Pull to Lock	103WL2	-	No Charge
<input type="checkbox"/> Pull to Lock + Wheel Lock Extension	103WL3	E0961	\$60.00


21. ANTI TIP			
<input type="checkbox"/> Quick FlipUp Anti Tipper	103MA23	E0971	\$120.00

22. POSITIONING STRAP			
<input type="checkbox"/> Auto Seat Positioning Strap	103MA03	E0978	\$52.00
<input type="checkbox"/> Velcro Release Positioning Strap	103MA05	E0978	\$48.00

23. MISCELLANEOUS ACCESSORIES			
<input type="checkbox"/> Quick Release Axle	103AL1	-	No Charge
<input type="checkbox"/> Multiback Cane Mount Lateral 4"x 6"	CML4	-	\$100.00
<input type="checkbox"/> Multiback Cane Mount Lateral 5"x 9"	CML5	-	\$100.00
<input type="checkbox"/> Multiback Cane Mount Lateral 6"x 10"	CML6	-	\$100.00
<input type="checkbox"/> Padded Lap Tray 15-16"	103MA09	E0950	\$230.00
<input type="checkbox"/> Padded Lap Tray 17-18"	103MA10	E0950	\$230.00
<input type="checkbox"/> Padded Lap Tray 19-20"	103MA11	E0950	\$230.00
<input type="checkbox"/> Half Padded Lap Tray, left side	103MA14	E0950	\$350.00
<input type="checkbox"/> Half Padded Lap Tray, right side	103MA15	E0950	\$350.00
<input type="checkbox"/> IV Pole, telescopic	105MA20	K0105	\$175.00



STF	REAR	CASTER	FORK	FORK	AXLE PLATE	STF	STEM	SIDE
FRONT	WHEEL	SIZE	SIZE	HOLE	HOLE	rear	SIZE	FRAME
(inches)	(inches)	(inches)				(inches)		STYLE
13.00	20	4	SHORT	A	SUPER LOW	12.50	SHORT	LOW
14.00	20	4	LRG	A	1	13.50	SHORT	LOW
14.50	20	4	LRG	B	1	14.00	SHORT	LOW
15.00	20	4	LRG	C	2	14.25	SHORT	LOW
15.00	20	5	LRG	B	2	14.25	SHORT	LOW
15.50	20	5	LRG	C	2	14.75	SHORT	LOW
15.50	22	4	LRG	B	2	14.75	LONG	LOW
16.00	20	4	LRG	C	3	15.25	LONG	LOW
16.00	22	5	LRG	B	2	15.25	LONG	LOW
16.50	20	5	LRG	C	3	16.00	LONG	LOW
17.00	20	6	LRG	B	3	16.25	SHORT	STANDARD
17.00	22	5	LRG	C	2	15.75	SHORT	STANDARD
17.50	20	6	LRG	C	4	17.00	SHORT	STANDARD
17.50	22	5	LRG	B	3	17.00	LONG	STANDARD
18.00	20	6	LRG	B	4	17.25	LONG	STANDARD
18.00	22	6	LRG	B	3	17.00	LONG	STANDARD
18.00	24	6	LRG	B	3	17.50	LONG	STANDARD
18.50	22	6	LRG	C	4	18.25	LONG	STANDARD
18.50	24	6	LRG	C	3	18.00	LONG	STANDARD
19.00	22	7	LRG	C	4	18.25	LONG	STANDARD
19.50	22	7	LRG	C	4	18.25		STANDARD



**WARRANTY**

Future Mobility Healthcare warrants the Galaxy Lite wheelchair to be free from defects in materials and workmanship for a period of one (1) year from the date of purchase. The frame is warranted for the lifetime of the original purchaser and/or user.