



FUTURE

MOBILITY HEALTHCARE INC.

ADDUCTOR PAD

PARTS MANUAL

FUTURE MOBILITY HEALTHCARE APPRECIATES YOUR BUSINESS.
OUR CUSTOMER SERVICE REPRESENTATIVES CAN BE CONTACTED AT:

**TO PLACE ORDERS:
(905) 671-1661**

FOR YOUR CONVENIENCE YOU MAY ALSO PLACE ORDERS BY FAX AT:

FAX: (905) 671-3377

PLEASE VISIT OUR WEBSITE AT:

www.futuremobility.ca

*Part prices and configurations included in this Parts Manual are subject to change without prior notice.

* All prices are shown as per unit individual price.



Contents

ADDUCTOR PAD MOUNTING HARDWARE 4


ADDUCTOR PAD MOUNTING HARDWARE ASSEMBLY..... 4

ADDUCTOR PAD MOUNTING HARDWARE

ADDUCTOR PAD MOUNTING HARDWARE ASSEMBLY



A) INSET MOUNT



B) FLUSH MOUNT

ADDUCTOR PAD MOUNTING HARDWARE ASSEMBLY - TABLE

ITEM	QTY	NOTES	PART NO.	DESCRIPTION	PRICE/ PER UNIT
1	4		1000397	NUT, HEX-NYLOCK THIN, 1/4-20NC	\$2.00
2	2		1000411	SCREW, SOC BUTTON HEAD 1/4-20 X 0.500	\$3.00
3	4		1001391	SCREW, SOCKET HEAD BUTTON, 1/4-20NC X 3/4	\$3.00
4	4		1004069	SCREW, SOC BUTTON HEAD 1/4-20 X 0.375	\$3.00
5	1	B	1021664	ADDUCTOR PAD EXTRA SMALL 5" WIDE W/PEM NUTS	\$30.00
	1	B	1013238	ADDUCTOR PAD SMALL 10" WIDE W/PEM NUTS	\$40.00
	1	B	1013239	ADDUCTOR PAD MEDIUM 12" WIDE W/PEM NUTS	\$40.00
	1	B	1013240	ADDUCTOR PAD LARGE 14" WIDE W/PEM NUTS	\$40.00
6	2		1013245	ADDUCTOR MOUNT BRACKET BTM	\$30.00
7	2		1013246	ADDUCTOR MOUNT BRKT FLUSH AT END	\$30.00
	2		1013247	ADDUCTOR MOUNT BRACKET OFFSET 1.25" AT END	\$30.00
8	4		1020550	1/4 X5/8 X 0.05 S.STEEL WASHER	\$1.50
	1	A	1038749	KIT, COMPLETE H'WARE AND MOUNTING BKTS FLUSH MOUNT EXTRA SMALL	\$185.00
	1	A	1038742	KIT, COMPLETE H'WARE AND MOUNTING BKTS FLUSH MOUNT SMALL	\$195.00
	1	A	1038743	KIT, COMPLETE H'WARE AND MOUNTING BKTS FLUSH MOUNT MEDIUM	\$195.00
	1	A	1038744	KIT, COMPLETE H'WARE AND MOUNTING BKTS FLUSH MOUNT LARGE	\$195.00
	1	A	1038748	KIT, COMPLETE H'WARE AND MOUNTING BKTS OFFSET MOUNT EXTRA SMALL	\$185.00
	1	A	1038745	KIT, COMPLETE H'WARE AND MOUNTING BKTS OFFSET MOUNT SMALL	\$195.00
	1	A	1038746	KIT, COMPLETE H'WARE AND MOUNTING BKTS OFFSET MOUNT MEDIUM	\$195.00
	1	A	1038747	KIT, COMPLETE H'WARE AND MOUNTING BKTS OFFSET MOUNT LARGE	\$195.00

NOTES:

A- EACH KIT COMES AS A SET OF 2, FOR BOTH SIDES OF THE CHAIR.

B- FOAM AND COVER NOT SHOWN