

LUNAR HEAVY DUTY™

The **Lunar Back Heavy Duty™** was designed for users weighing up to 500lbs who require aggressive lateral support and more complex positioning. Reinforcement ribs and extra-large shell allow this back to maintain its function in any high load condition. Lunar's sleek and attractive tailored look will appeal to a wide variety of users!



Lunar Heavy Duty

STYLE: Mid Support

Stock Sizes – Same Day Shipment

Part Number	Size W x H	Description
MBSHDBL2216	22" x 16"	Basic
MBSHDUL2216	22" x 16"	Ultra
MBSHDBL2218	22" x 18"	Basic
MBSHDUL2218	22" x 18"	Ultra
MBSHDBL2416	24" x 16"	Basic
MBSHDUL2416	24" x 16"	Ultra
MBSHDBL2418	24" x 18"	Basic
MBSHDUL2418	24" x 18"	Ultra

Part Number	Size W x H	Description
MBSHDBL2616	26" x 16"	Basic
MBSHDUL2616	26" x 16"	Ultra
MBSHDBL2618	26" x 18"	Basic
MBSHDUL2618	26" x 18"	Ultra
MBSHDBL2816	28" x 16"	Basic
MBSHDUL2816	28" x 16"	Ultra
MBSHDBL2818	28" x 18"	Basic
MBSHDUL2818	28" x 18"	Ultra

Part Number	Size W x H	Description
MBSHDBL3016	30" x 16"	Basic
MBSHDUL3016	30" x 16"	Ultra
MBSHDBL3018	30" x 18"	Basic
MBSHDUL3018	30" x 18"	Ultra
MBSHDBL3216	32" x 16"	Basic
MBSHDUL3216	32" x 16"	Ultra
MBSHDBL3218	32" x 18"	Basic
MBSHDUL 3218	32" x 18"	Ultra

Lunar Heavy Duty

STYLE: Full Support

Stock Sizes – Same Day Shipment

Part Number	Size W x H	Description
FBSHDBL2220	22" x 20"	Basic
FBSHDUL2220	22" x 20"	Ultra
FBSHDBL2222	22" x 22"	Basic
FBSHDUL2222	22" x 22"	Ultra
FBSHDBL2224	22" x 24"	Basic
FBSHDUL2224	22" x 24"	Ultra
FBSHDBL2420	24" x 20"	Basic
FBSHDUL2420	24" x 20"	Ultra
FBSHDBL2422	24" x 22"	Basic
FBSHDUL 2422	24" x 22"	Ultra
FBSHDBL2424	24" x 24"	Basic
FBSHDUL2424	24" x 24"	Ultra

Part Number	Size W x H	Description
FBSHDBL2620	26" x 20"	Basic
FBSHDUL2620	26" x 20"	Ultra
FBSHDBL2622	26" x 22"	Basic
FBSHDUL2622	26" x 22"	Ultra
FBSHDBL2624	26" x 24"	Basic
FBSHDUL2624	26" x 24"	Ultra
FBSHDBL2820	28" x 20"	Basic
FBSHDUL2820	28" x 20"	Ultra
FBSHDBL2822	28" x 22"	Basic
FBSHDUL2822	28" x 22"	Ultra
FBSHDBL2824	28" x 24"	Basic
FBSHDUL2824	28" x 24"	Ultra

Part Number	Size W x H	Description
FBSHDBL3020	30" x 20"	Basic
FBSHDUL3020	30" x 20"	Ultra
FBSHDBL3022	30" x 22"	Basic
FBSHDUL3022	30" x 22"	Ultra
FBSHDBL3024	30" x 24"	Basic
FBSHDUL3024	30" x 24"	Ultra
FBSHDBL3220	32" x 20"	Basic
FBSHDUL3220	32" x 20"	Ultra
FBSHDBL3222	32" x 22"	Basic
FBSHDUL3222	32" x 22"	Ultra
FBSHDBL3224	32" x 24"	Basic
FBSHDUL3224	32" X 24"	Ultra

PRODUCT HIGHLIGHTS

- Lightweight, sturdy aluminum shell
- Reinforcement ribs
- Attractive cover design
- Replaceable foam
- Available in 2 lateral styles (basic lateral and ultra lateral)
- Available in 2 back style heights (mid-back and full back)
- Optional lateral supports, available in 4 styles: Fixed, Height Adjustable, Fixed Offset and Adjustable Swing Away
- Pins adjustable to accommodate up to 2" difference in seat width
- Fully adjustable for height, depth and angle of recline
- Hardware fits 3/4", 7/8" and 1" tubing and can be installed within a few minutes with the right tools
- Quick release lower latch levers allows the backrest to be easily removed.
- **Product Weight:** 8.55 lbs (22"W x 18"H Basic)
Capacity: 500 lbs
Warranty: 3 years limited on back and 90 days on cover, 1 year on hardware



Readily available for headrest





FUTURE MOBILITY

CANADA

6750 Professional Court, Mississauga,
Ontario, L4V 1X6



1-888-737-4011



www.futuremobility.ca



USA



11126 Greenstone Ave, Santa Fe Springs
CA 90670



716-783-9130



www.futuremobility.com

