



AIR BACK



Lightweight aluminum shell

6" of lateral support



Raised lateral for high trunk support

Air bladder helps with lateral positioning provides further lumbar support

Max Series Backs **STYLE: Max Air Back** Stock Sizes – Same Day Shipment

Part Number	Size W x H	Part Number	Size W x H
MAB1616	16" x 16"	MAB1820	18" x 20"
MAB1620	16" x 20"	MAB2016	20" x 16"
MAB1816	18" x 16"	MAB2020	20" x 20"

Non-Stock Sizes – Shipment within 3 business days

Part Number	Size W x H	Part Number	Size W x H
MAB1414	14" x 14"	MAB1818	18" x 18"
MAB1618	16" x 18"	MAB1916	19" x 16"
MAB1716	17" x 16"	MAB1920	19" x 20"
MAB1720	17" x 20"	MAB2018	20" x 18"

Custom Sizes - Shipment within 7 business days

Part Number	Size W x H
MABCUS	Custom

Options

Part Number	Description
MABLT	Omit Lateral Inserts
MABLB	Omit Lumbar Inserts

Bariatric Sizes – Shipment within 3 business days

Part Number	Size W x H	Part Number	Size W x H
MAB2116	21" x 16"	MAB2316	23" x 16"
MAB2118	21" x 18"	MAB2318	23" x 18"
MAB2120	21" x 20"	MAB2320	23" x 20"
MAB2216	22" x 16"	MAB2416	24" x 16"
MAB2218	22" x 18"	MAB2418	24" x 18"
MAB2220	22" x 20"	MAB2420	24" x 20"


Custom Bariatric Sizes - Shipment within 7 business days

Part Number	Size W x H
MABBCUS	Custom

PRODUCT HIGHLIGHTS


- Fully adjustable for height, depth and angle of recline
- Hardware fits 3/4", 7/8" and 1" tubing and can be installed within a few minutes with the right tools
- Quick release lower latch levers allows the backrest to be easily removed.
- Product Weight:** 5.7 lbs (18"W x 18"H)
Warranty: 2 years limited on back and 90 days on cover, 1 year on hardware

Lateral Dimensions		Max Air Back
MAX AIR, Heights 14, 16, 18, 20		
Chair Width	DIM B (")	DIM A (")
14"	10.5	7.6
16"	12.5	7.6
18"	14.5	7.6
20"	16.5	7.6
22"	18.5	7.6



Cover and Foam Layers, showing depth of the back

Standard Mounting Hardware





FUTURE MOBILITY

CANADA

6750 Professional Court, Mississauga,
Ontario, L4V 1X6



1-888-737-4011



www.futuremobility.ca



USA



11126 Greenstone Ave, Santa Fe Springs
CA 90670



716-783-9130



www.futuremobility.com

