



TRUEFITT BACK



Available in Visco Foam
and Poly Foam

Max Series Backs **STYLE: Max TrueFitt Back**
Stock Sizes – Same Day Shipment

Part Number	Size W x H	Part Number	Size W x H
MTFB1616	16" x 16"	MTFB1820	18" x 20"
MTFB1620	16" x 20"	MTFB2016	20" x 16"
MTFB1816	18" x 16"	MTFB2020	20" x 20"

Non-Stock Sizes – Shipment within 3 business days

Part Number	Size W x H	Part Number	Size W x H
MTFB1414	14" x 14"	MTFB1818	18" x 18"
MTFB1618	16" x 18"	MTFB1916	19" x 16"
MTFB1716	17" x 16"	MTFB1920	19" x 20"
MTFB1720	17" x 20"	MTFB2018	20" x 18"

Custom Sizes - Shipment within 7 business days

Part Number	Size W x H
MTFBCUS	Custom

Options		
Part Number	Size W x H	Description
MTFBV	2" x 2"	Visco Foam
MTFB3	3" x 2"	Poly Foam

PRODUCT HIGHLIGHTS



- Fully adjustable for height, depth and angle of recline
- Hardware fits 3/4", 7/8" and 1" tubing and can be installed within a few minutes with the right tools
- Quick release lower latch levers allows the backrest to be easily removed.
- Product Weight:** 5.5 lbs (18"W x 16"H)
Warranty: 2 years limited on back and 90 days on cover, 1 year on hardware

Max Series Backs **STYLE: Max TrueFitt Back**
Bariatric Sizes – Shipment within 3 - 7 business days


Part Number	Size W x H	Part Number	Size W x H
MTFB2116	21" x 16"	MTFB2316	23" x 16"
MTFB2118	21" x 18"	MTFB2318	23" x 18"
MTFB2120	21" x 20"	MTFB2320	23" x 20"
MTFB2216	22" x 16"	MTFB2416	24" x 16"
MTFB2218	22" x 18"	MTFB2418	24" x 18"
MTFB2220	22" x 20"	MTFB2420	24" x 20"

Custom Bariatric Sizes - Shipment within 3 - 10 business days

Part Number	Size W x H
MTFBBCUS	Custom



Lateral Dimensions		Max Truefitt
MAX TRUEFITT, Heights 14, 16, 18, 20		
Chair Width	DIM B (")	DIM A (")
14"	10.5	7.6
16"	12.5	7.6
18"	14.5	7.6
20"	16.5	7.6
22"	18.5	7.6



2



1



3



4

Max TrueFitt Back showing:
1. side, 2. top, 3. bottom, 4. backside
of the foam layers and hardware



FUTURE MOBILITY

CANADA

6750 Professional Court, Mississauga,
Ontario, L4V 1X6



1-888-737-4011



www.futuremobility.ca



USA



11126 Greenstone Ave, Santa Fe Springs
CA 90670



716-783-9130



www.futuremobility.com

