



PRISM ADJUSTABLE TENSION BACK



PRODUCT HIGHLIGHTS

- Attractive design
- Pocket for personal items
- Two fixed straps and two removable straps
- All day comfort and pressure relief
- **Product Weight:** 2.6 lbs (18"W x 18"H)
Capacity: 250 lbs
Warranty: 1 years limited

The **Adjustable Tension Back™** is made for additional comfort and adjustability. Two fixed straps and two Velcro straps allow the user to apply tension where it's needed the most. Allows for all-day comfort and pressure relief for the wheelchair's end user.

Prism Series Backs STYLE: Adjustable Tension Back Tension Back Low - Same Day Shipment

Part Number	Size W x H	Part Number	Size W x H
TB1413	14" x 13"	TB1813	18" x 13"
TB1513	15" x 13"	TB1913	19" x 13"
TB1613	16" x 13"	TB2013	20" x 13"
TB1713	17" x 13"		

Tension Back High - Shipment within 3 business days

Part Number	Size W x H	Part Number	Size W x H
TB1416	14" x 16"	TB1816	18" x 16"
TB1516	15" x 16"	TB1916	19" x 16"
TB1616	16" x 16"	TB2016	20" x 16"
TB1716	17" x 16"		

Custom Sizes - Shipment within 7 business days

Part Number	Size W x H
TBCUS	Custom

Nylon upholstery that

4 individual adjustable straps for

Pocket for stor-





FUTURE MOBILITY

CANADA

6750 Professional Court, Mississauga,
Ontario, L4V 1X6



1-888-737-4011



www.futuremobility.ca



USA



11126 Greenstone Ave, Santa Fe Springs
CA 90670



716-783-9130



www.futuremobility.com

