



# PRISM BASIC BACK

The **Prism Basic Back™** was designed for users who require mild to moderate postural support. Extras, such as the positioning options, the ease of attachment and the unique method of adjusting for a wide range of angles, makes the Basic Back appealing for a wide variety of users.



### Prism Series Backs      STYLE: Basic Back Stock Sizes – Same Day Shipment

Part Number	Size W x H	Part Number	Size W x H
B1616	16" x 16"	B1820	18" x 20"
B1620	16" x 20"	B2016	20" x 16"
B1816	18" x 16"	B2020	20" x 20"

### Non-Stock Sizes – Shipment within 3 business days

Part Number	Size W x H	Part Number	Size W x H
B1414	14" x 14"	B1818	18" x 18"
B1618	16" x 18"	B1916	19" x 16"
B1716	17" x 16"	B1920	19" x 20"
B1720	17" x 20"		

### Custom Sizes - Shipment within 7 business days

Part Number	Size W x H
BCUS	Custom

### Prism Series Backs      STYLE: Basic Back Bariatric Sizes – Shipment within 3 business days

Part Number	Size W x H	Part Number	Size W x H
B2116	21" x 16"	B2316	23" x 16"
B2118	21" x 18"	B2318	23" x 18"
B2120	21" x 20"	B2320	23" x 20"
B2216	22" x 16"	B2416	24" x 16"
B2218	22" x 18"	B2418	24" x 18"
B2220	22" x 20"	B2420	24" x 20"

### Custom Bariatric Sizes - Shipment within 7 business days


Part Number	Size W x H
BBCUS	Custom

# ▶ PRODUCT HIGHLIGHTS

- Lightweight, sturdy 5052 grade aluminum shell
- A full 5" range of forward and rear adjustable
- Back angle adjustment forward 15 degrees, back 20 degrees, optimizing the users' positioning
- Attractive cover design
- Pins adjustable to accommodate up to 2" difference in seat width (standard, 0.5" and 1" longer)
- Cane clamps hardware included
- Replaceable foam
- Quick latch release for easy installation
- Easily adjusted while the client stays in the chair
- **Product Weight:** 4.75 lbs (6.7 lbs bariatric), (18"W x 16"H)  
**Capacity:** 250 lbs (350 lbs bariatric)  
**Warranty:** 3 years limited on back and 90 days on cover, 1 year on hardware



## Standard Mounting Hardware





# FUTURE MOBILITY

## CANADA

6750 Professional Court, Mississauga,  
Ontario, L4V 1X6



1-888-737-4011



[www.futuremobility.ca](http://www.futuremobility.ca)



## USA



11126 Greenstone Ave, Santa Fe Springs  
CA 90670



716-783-9130



[www.futuremobility.com](http://www.futuremobility.com)

