



PRISM SUPPER BACK FOR KIDS

The **Super Back for Kids™** is a Junior size Back created with three layers of foam to give children the best in pressure relief. An anatomically designed foam base in this model gives superior back stability and positioning for the end user.



Prism Series Backs **STYLE: Super Back for Kids**
Stock Sizes - Same Day Shipment

Part Number	Size W x H	Part Number	Size W x H
KSB1012	10" x 12"	KSB1414	14" x 14"
KSB1212	12" x 12"	KSB1616	16" x 16"

Non-Stock Sizes – Shipment within 3 business days

Part Number	Size W x H	Part Number	Size W x H
KSB1112	11" x 12"	KSB1412	14" x 12"
KSB1114	11" x 14"	KSB1416	14" x 16"
KSB1214	12" x 14"	KSB1512	15" x 12"
KSB1216	12" x 16"	KSB1514	15" x 14"
KSB1312	13" x 12"	KSB1612	16" x 12"
KSB1316	13" x 16"	KSB1614	16" x 14"

Custom Sizes - Shipment within 7 business days

Part Number	Size W x H
KSBCUS	Custom

Back view of Foam Layer and Aluminum Shell

PRODUCT HIGHLIGHTS

- Lightweight, sturdy aluminum 5052 grade shell
- Full 5" range of forward and rear adjustment (angle adjustment forward 15° and back 20°)
- Pins adjustable to accommodate up to 2" difference in seat width (standard, 0.5" and 1" longer)
- Attractive cover design
- Cane clamps hardware included
- Replaceable foam
- Lumbar pad and contour pad are standard
- Quick latch release for easy installation
- Fully adjustable for height, depth and angle of recline
- Hardware fits 3/4", 7/8" and 1" tubing and can be installed within a few minutes with the right tools
- Quick release lower latch levers allows the backrest to be easily removed.
- **Product Weight:** 3.2 lbs (12"W x 14"H)
Capacity: 150 lbs
Warranty: 3 years limited on back and 90 days on cover, 1 year on hardware




Cover, Foam Layer and Aluminum Hardware

Standard Mounting Hardware



Lateral Dimensions		Kids Back
KIDS BACK, Heights 10, 12, 13, 14, 16		
Chair Width	DIM B (")	DIM A (")
12"	9.6	4.5
13"	10.6	4.5
14"	11.6	4.5
15"	12.6	4.5
16"	13.6	4.5





FUTURE MOBILITY

CANADA

6750 Professional Court, Mississauga,
Ontario, L4V 1X6



1-888-737-4011



www.futuremobility.ca



USA



11126 Greenstone Ave, Santa Fe Springs
CA 90670



716-783-9130



www.futuremobility.com

