



FUTURE

MOBILITY PRODUCTS INC.

DIVISION OF ABF MFG GROUP INC.
WWW.FUTUREMOBILITY.COM



TEL: 888-342-3119
FAX: 424-217-7327
orderdesk@futuremobility.com

Account Information

Account #: _____
Business Name: _____
Ship to Address: _____
Contact: _____
Tag For: _____

Order Date: _____
P.O.# _____
Phone: _____
City: _____
State: _____ Zip Code: _____

The HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions. Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy. Prices are subject to change without notice

1. FRAME TYPE

	Parts #	E-Code	Retail Price
<input type="checkbox"/> SuperChair For Kids	104	E1233	\$2,765.00

2. FRAME FINISH

<input type="checkbox"/> Wet Black		104C01	-	No Charge
<input type="checkbox"/> Silver		104C02	-	No Charge
<input type="checkbox"/> Blue		104C03	-	No Charge
<input type="checkbox"/> Green		104C04	-	No Charge
<input type="checkbox"/> Orange		104C05	-	No Charge
<input type="checkbox"/> Purple		104C06	-	No Charge
<input type="checkbox"/> Red		104C07	-	No Charge
<input type="checkbox"/> Yellow		104C08	-	No Charge
<input type="checkbox"/> Custom Colors		104C17	-	\$300.00

3. SEAT WIDTH

<input type="checkbox"/> 10"	104SW10	-	No Charge
<input type="checkbox"/> 11"	104SW11	-	No Charge
<input type="checkbox"/> 12"	104SW12	-	No Charge
<input type="checkbox"/> 13"	104SW13	-	No Charge
<input type="checkbox"/> 14"	104SW14	-	No Charge
<input type="checkbox"/> 15"	104SW15	-	No Charge
<input type="checkbox"/> 16"	104SW16	-	No Charge

4. SEAT DEPTH

<input type="checkbox"/> 11"	104SD11	-	No Charge
<input type="checkbox"/> 12"	104SD12	-	No Charge
<input type="checkbox"/> 13"	104SD13	-	No Charge
<input type="checkbox"/> 14"	104SD14	-	No Charge
<input type="checkbox"/> 15"	104SD15	-	No Charge
<input type="checkbox"/> 16"	104SD16	-	No Charge

5. SEAT TO FLOOR HEIGHT

<input type="checkbox"/> 17.5" Seat to Floor (20" Wheels & 4" Casters)	104STF62	-	No Charge
<input type="checkbox"/> 18" Seat to Floor (20" Wheels & 4" Casters)	104STF63	-	No Charge
<input type="checkbox"/> 18.5" Seat to Floor (22" Wheels & 5" Casters)	104STF64	-	No Charge
<input type="checkbox"/> 19" Seat to Floor (22" Wheels & 7" Casters)	104STF60	-	No Charge
<input type="checkbox"/> 19.5" Seat to Floor (22" Wheels & 7" Casters)	104STF61	-	No Charge
<input type="checkbox"/> 19.75" Seat to Floor (22" Wheels & 8" Casters)	104STF65	-	No Charge
<input type="checkbox"/> 18.5" Seat to Floor (12" Wheels & 5" Casters)	104STF66	-	No Charge
<input type="checkbox"/> 19" Seat to Floor (12" Wheels & 7" Casters)	104STF67	-	No Charge
<input type="checkbox"/> 19.5" Seat to Floor (12" Wheels & 7" Casters)	104STF68	-	No Charge
<input type="checkbox"/> 20" Seat to Floor (12" Wheels & 8" Casters)	104STF69	-	No Charge

6. REAR WHEEL RIM

<input type="checkbox"/> Mag Wheels	104RW09	-	No Charge
<input type="checkbox"/> Spoke Wheels	104RW10	-	\$75.00

7. REAR WHEEL TIRE

<input type="checkbox"/> Airless (Solid Rubber)	104RT1	-	No Charge
<input type="checkbox"/> Pneumatic (Only Available for 20" or 22")	104RT2	E2211	\$90.00

SUPERKIDS FRONT SEAT TO FLOOR HEIGHT MATRIX

Front STF	Rear Wheel	Caster Size	Fork	Fork Hole	Axel Hole	Rear STF
18.5"	12"	4"	Long	C	A	17.25"
18.5"	12"	5"	Long	B	A	17.25"
19"	12"	7"	Short	C	A	17.5"
19.5"	12"	7"	Long	B	B	18.25"
20"	12"	8"	Long	B	B	18.5"
17.5"	20"	4"	Short	C	A	15.75"
18"	20"	5"	Short	C	A	15.75"
18.5"	22"	5"	Long	B	A	16.75"
19"	20"	7"	Short	C	B	17"
19.5"	22"	7"	Long	B	B	18"
20"	22"	7"	Long	C	B	18"

SuperKids Front STF Height Availability							
Rear Wheel	4" Caster						
20"	17.5"	18"					
Rear Wheel	5" Caster						
12"			18.5"				
22"			18.5"				
Rear Wheel	7" Caster						
12"			19"	19.5"			
22"			19"	19.5"			
Rear Wheel	8" Caster						
12"					20"		
22"					19.75"		


8. HANDRIM

	Parts #	E-Code	Retail Price
<input type="checkbox"/> Chrome	104HR1	-	No Charge
<input type="checkbox"/> Plastic Coated	104HR2	K0059	\$160.00
<input type="checkbox"/> No Handrims	104HR3	-	No Charge

9. ARM TYPE

<input type="checkbox"/> T-Style Adjustable Height Pin Release	104AR05	E0973	No Charge
<input type="checkbox"/> T-Style Fixed Height	104AR02	E0973	No Charge
<input type="checkbox"/> Cantilever Arm Adj. Angle	104AR11	E0973	\$225.00
<input type="checkbox"/> Cantilever Arm Adj. Angle with Guard	104AR12	E0973	\$320.00

Armrest Matrix for T - Style Pin Release Arm						
Armrest Style	Transfer Tube Length D1	Arm Guard Length D2	H Range	Arm Pad Length L		
				Desk	Full	Curved
Youth	4.50"	5.00"	8.75" - 13.75"	10.25"	14.0"	17.0"



10. ARM PAD TYPE

	Parts #	E-Code	Retail Price
<input type="checkbox"/> Desk length Vinyl 10" (pair)	104AP1	-	No Charge
<input type="checkbox"/> Full length Vinyl 14" (pair)	104AP2	-	No Charge
<input type="checkbox"/> Contour Gel 17" (pair)	104AP3	-	No Charge

11. BACK POST

<input type="checkbox"/> Angle Adj 10° Bend w/Push Handle (low)	104SP3	-	No Charge
<input type="checkbox"/> Angle Adj 10° Bend w/Push Handle (high)	104SP4	-	\$90.00
<input type="checkbox"/> Angle Adj 15° Bend w/Push Handle (low)	104PC3	-	No Charge
<input type="checkbox"/> Angle Adj 15° Bend w/Push Handle (high)	104PC4	-	\$90.00

12. BACK POST HEIGHT

<input type="checkbox"/> 10"	NA	-	No Charge
<input type="checkbox"/> 12"	NA	-	No Charge
<input type="checkbox"/> 14"	NA	-	No Charge
<input type="checkbox"/> 16"	NA	-	No Charge

13. BACK OPTIONS

<input type="checkbox"/> No Back	104BK49	-	Standard
<input type="checkbox"/> Super Chair Kids Back Stock (10",12",14",16")	104BK75	E2620	\$695.00
<input type="checkbox"/> Super Chair Kids Back Non Stock (11" - 16")	104BK76	E2620	\$695.00
<input type="checkbox"/> Super Chair Kids Back Custom	104BK77	E2620	\$745.00
<input type="checkbox"/> Slip-on Padded Back Upholstery	104BK48	-	\$150.00
<input type="checkbox"/> Prism Basic Back 14" - 20" wide	BB14-20	E2620	\$615.00
<input type="checkbox"/> Prism Ultra Back 14" - 20" wide	UB14-20	E2620	\$615.00
<input type="checkbox"/> Prism Air Back 14" - 20" wide	AB14-20	E2615	\$655.00
<input type="checkbox"/> Prism TrueFitt Back 14" - 20" wide	TFB14-20	E2615	\$625.00
<input type="checkbox"/> Max Ultra Back 14" - 20" wide	MUB14-20	E2620	\$645.00
<input type="checkbox"/> Max Air Back 14" - 20" wide	MAB14-20	E2620	\$725.00
<input type="checkbox"/> Max TrueFitt Back 14" - 20" wide	MTFB14-20	E2620	\$675.00

W () x H ()

14. BACK PINS

<input type="checkbox"/> Not Required	NA	-	No Charge
<input type="checkbox"/> Longer Pins for Back Hardware 1" Growth	104BP1	-	\$50.00
<input type="checkbox"/> Longer Pins for Back Hardware 2" Growth	104BP2	-	\$50.00

15. SEAT PAN OPTIONS

<input type="checkbox"/> Not Required	NA	-	No Charge
<input type="checkbox"/> Expandable Seat Pan	104SE01	E2231	\$150.00
<input type="checkbox"/> Prism Drop Base FM	104SC36	E2231	\$100.00
<input type="checkbox"/> Prism Drop Base 1"	104SC37	E2231	\$100.00
<input type="checkbox"/> Prism Drop Base 2"	104SC38	E2231	\$100.00
<input type="checkbox"/> Prism Adjustable Drop Base	104SC42	E2231	\$150.00

16. CUSHION OPTIONS

<input type="checkbox"/> No Cushion	104SC20	-	No Charge
<input type="checkbox"/> Super Kids Cushion Stock	104SC44	E2231	\$150.00
<input type="checkbox"/> Super Kids Cushion Non Stock	104SC45	E2601	\$100.00
<input type="checkbox"/> Super Kids Cushion Custom	104SC46	E2605	\$300.00
<input type="checkbox"/> Prism Ideal	104SC03	E2601	\$100.00
<input type="checkbox"/> Prism Elite	104SC06	E2601	\$175.00
<input type="checkbox"/> Prism Supreme	104SC09	E2603	\$225.00
<input type="checkbox"/> Prism Supreme Air	104SC15	E2603	\$300.00
<input type="checkbox"/> Prism Supreme II Plus	104SC12	E2607	\$365.00

W () x D ()

17. FRONT RIGGING

<input type="checkbox"/> 70° Swingaway	104FR03	-	No Charge
<input type="checkbox"/> 80° Swingaway	104FR19	-	\$97.00
<input type="checkbox"/> 90° Swingaway	104FR05	-	\$97.00
<input type="checkbox"/> Swingaway Elevating Legrest (pair)	104FR10	E0990	\$275.00
<input type="checkbox"/> No Front Rigging	104FR11	-	No Charge

18. FOOTPLATE

<input type="checkbox"/> Super Chair Aluminum Flip-Up - Small	104FP13	-	No Charge
<input type="checkbox"/> Super Chair Aluminum Flip-Up - Medium	104FP14	-	No Charge
<input type="checkbox"/> Super Chair Aluminum Flip-Up - Large	104FP15	-	No Charge
<input type="checkbox"/> Composite Flip Up - Small	104FP01	-	No Charge
<input type="checkbox"/> Aluminum Flip Up - Medium	104FP06	-	\$48.00
<input type="checkbox"/> Aluminum Angle Adj Flip Up Footplates - Small	104FP10	K0040	\$120.00

WARRANTY:

Future Mobility Healthcare warrants the e SuperChair for Kids to be free from defects in materials and workmanship for a period of one (1) year from the date of purchase. The frame is warranted for the lifetime of the original purchaser and/or user.

19. FOOTPLATE / LEGREST ACCESSORIES

	Parts #	E-Code	Retail Price
<input type="checkbox"/> Heel Loops	104FA01	K0951	\$40.00
<input type="checkbox"/> No Heel Loops	104FA02	-	No Charge
<input type="checkbox"/> Padded Calf Strap	104FA03	K0038	\$45.00
<input type="checkbox"/> Foam Pad on Footplate (pair)	104FA07	-	\$45.00

20. WHEEL LOCKS

<input type="checkbox"/> Push to Lock	104WL1	-	No Charge
<input type="checkbox"/> Push to Lock + Wheel Lock Extension	104WL3	E0961	\$60.00
<input type="checkbox"/> Pull to Lock	104WL2	-	No Charge
<input type="checkbox"/> Pull to Lock + Wheel Lock Extension	104WL3	E0961	\$60.00

21. HEADREST HARDWARE

<input type="checkbox"/> Omit Headrest Hardware	104HE99	-	Standard
<input type="checkbox"/> Fixed Hardware	104HE01	E0955	\$195.00
<input type="checkbox"/> Fixed Hardware (+ 4")	104HE02	E0955	\$200.00
<input type="checkbox"/> Single Multi Axial	104HE03	E0955	\$217.00
<input type="checkbox"/> Single Multi Axial (+ 4")	104HE04	E0955	\$222.00
<input type="checkbox"/> Dual Multi Axial	104HE05	E0955	\$219.00
<input type="checkbox"/> Dual Multi Axial (+ 4")	104HE06	E0955	\$225.00
<input type="checkbox"/> Links Multi Adj. Hardware (Comes with 2 links)	104HE10	E0955	\$240.00
<input type="checkbox"/> Links Multi Adj. Curved Hardware (Comes with 2 links)	104HE11	E0955	\$240.00
<input type="checkbox"/> Additional Link for Headrest	104HE12	E0955	\$50.00
<input type="checkbox"/> 2 Additional Links for Headrest	104HE13	E0955	\$100.00
<input type="checkbox"/> 3 Additional Links for Headrest	104HE14	E0955	\$150.00

22. HEADREST ADAPTOR BRACKETS

<input type="checkbox"/> Universal Headrest Bracket	104HEBA	E1028	\$86.00
---	---------	-------	---------

23. HEADREST FOAM PAD

<input type="checkbox"/> Omit Headrest Foam Pad	104HE09	-	Standard
<input type="checkbox"/> Kids Headrest Foam Pad	104HF5	-	No Charge
<input type="checkbox"/> Small Headrest Foam Pad	104HF1	E0955	No Charge
<input type="checkbox"/> Small Headrest with Wings Foam Pad	104HF2	E0955	\$90.00
<input type="checkbox"/> Medium Headrest Foam Pad	104HF7	E0955	No Charge
<input type="checkbox"/> Medium Headrest with Wings Foam Pad	104HF8	E0955	\$90.00
<input type="checkbox"/> Large Headrest Foam Pad	104HF3	E0955	No Charge
<input type="checkbox"/> Large Headrest with Wings Foam Pad	104HF4	E0955	\$95.00
<input type="checkbox"/> Pro-Fit - Small	104HF12	E0955	\$175.00
<input type="checkbox"/> Pro-Fit - Medium	104HF13	E0955	\$195.00
<input type="checkbox"/> Pro-Fit - Large	104HF14	E0955	\$205.00

24. ANTI TIP

<input type="checkbox"/> Fixed	104MA22	E0971	\$90.00
<input type="checkbox"/> Quick FlipUp	104MA23	E0971	\$120.00

25. POSITIONING STRAP

<input type="checkbox"/> Auto	104MA03	E0978	\$52.00
<input type="checkbox"/> Airline	104MA04	E0978	\$76.00
<input type="checkbox"/> Velcro	104MA05	E0978	\$48.00

26. MISCELLANOUS ACCESSORIES

<input type="checkbox"/> Quick Release Axle	104AL1	-	\$60.00
<input type="checkbox"/> Wheelchair Attachment Tie Down	104MA58	-	\$150.00
<input type="checkbox"/> Foam Pad on Arm Guard	104MA06	E0978	\$76.00
<input type="checkbox"/> Padded Lap Tray 15-16"	104MA09	E0950	\$230.00
<input type="checkbox"/> Half Padded Lap Tray, left side	104MA14	E0950	\$350.00
<input type="checkbox"/> Half Padded Lap Tray, right side	104MA15	E0950	\$350.00
<input type="checkbox"/> IV Pole, Telescopic	104MA20	K0105	\$175.00
<input type="checkbox"/> Oxygen Tank Holder	104MA21	E2208	\$175.00
<input type="checkbox"/> Spoke Guards available on 20" or 22" wheels	104MA49	-	\$130.00

NOTES: